

User Guide for Scanalyze AI

Team 16, Ace
Harrison Chiu
Hamza Issa
Ahmad Hamadi
Jared Paul
Gurnoor Bal

April 4, 2025

Contents

| | |
|--|-----------|
| 1 Introduction | 5 |
| Overview of the Scanalyze Imaging Platform and its intended use by radiologists. | |
| 2 Signup | 6 |
| Instructions for creating a new user account, including required fields and system validation. | |
| 2.1 How to Sign Up | 6 |
| 2.2 What Happens Next | 6 |
| 3 Login | 7 |
| Steps to log in to the platform using your registered credentials. | |
| 3.1 How to Log In | 7 |
| 3.2 What Happens Next | 7 |
| 3.3 Important Notes | 8 |
| 4 Select Model | 8 |
| Guide to choosing the appropriate AI model (Chest X-Ray Analysis, Brain MRI Segmentation, Bone Fracture Detection) and proceeding to analysis. | |
| 4.1 Available Models | 8 |
| 4.2 How to Proceed | 8 |
| 5 Upload CXR Image | 10 |
| Instructions for uploading chest X-ray images using drag-and-drop or file selection. | |
| 5.1 How to Upload an Image | 10 |
| 5.2 Important Notes | 10 |
| 6 Predicted Results and Report of Findings | 12 |
| Detailed breakdown of the AI-generated diagnostic report | |
| 6.1 Report Overview | 12 |
| 6.2 Confidence Slider | 12 |
| 6.3 Findings | 12 |
| 6.4 Probability Distribution Table | 13 |
| 6.5 Heatmap Visualization | 13 |

7 Exploring the CXR Image 15

Interactive features for zooming and dragging the image to explore AI-highlighted findings.

- 7.1 Interactive Tools for Image Exploration 15
- 7.2 Practical Use 15

8 Download Report 17

How to generate and save a complete PDF report including findings, descriptions, probability scores, heatmap, and confidence settings.

- 8.1 Purpose of the Report 17
- 8.2 What the Report Includes 17
- 8.3 Using the Button 17

9 Frequently Asked Questions (FAQ) 18

Common questions regarding usage, data formats, diagnostic value, and technical features of the platform.

- 9.1 General Questions 18
 - 9.1.1 What is the purpose of the Scanalyze Imaging Platform? 18
 - 9.1.2 Is this platform a replacement for a radiologist's diagnosis? 18
- 9.2 Usage Questions 18
 - 9.2.1 What image formats are supported for upload? 18
 - 9.2.2 How long does it take to analyze an image? 18
 - 9.2.3 Can I use the platform on mobile devices? 18
- 9.3 Technical Questions 18
 - 9.3.1 What happens if the platform cannot process my image? 18
 - 9.3.2 How secure is my data? 19
 - 9.3.3 Can I download the analysis results? 19
- 9.4 AI and Diagnostic Questions 19
 - 9.4.1 How accurate are the AI predictions? 19
 - 9.4.2 What diseases can the AI detect? 19
 - 9.4.3 What should I do if the AI gives low-confidence predictions? 19

Revision History

| Date | Author | Section |
|---------|----------|--|
| 4 April | Hamza | All sections, finished creating User Guide |
| 4 April | Harrison | Edited all sections, added Login and FAQs |

1 Introduction

Welcome to the user guide for the Scanalyze Imaging Platform, an advanced AI-powered tool designed specifically for radiologists. This platform allows medical professionals to upload chest X-ray (CXR) images and receive automated diagnostic support, including disease prediction, comprehensive report generation, and interactive visual analysis tools.

The goal of this platform is to assist radiologists in their workflow by providing accurate, explainable outputs that complement professional clinical judgment. While this tool leverages advanced machine learning algorithms to generate insights, it is intended for use as a decision-support system, not a replacement for a formal radiologic diagnosis.


This guide covers the full range of platform functionality, from account creation and image upload to interpreting AI-generated reports, exploring image findings, and downloading clinical summaries.

2 Signup

The first step to using the platform is creating your account.

2.1 How to Sign Up

On the signup page, you will see a registration form like the one shown below:

A screenshot of a web form for 'Scanalyze Imaging'. The form has a light blue header with the logo and name. Below the header, there are three input fields: 'Username' with the value 'steven', 'Email' with the value 'steven@gmail.com', and 'Password' with masked characters '*****'. A blue 'Sign Up' button is positioned below the password field. At the bottom, there is a link that says 'Already have an account? Sign in'.

To complete the signup process:

1. Enter your Username – this will identify you on the platform.
2. Enter a valid Email address – make sure the email is correctly formatted (e.g., `name@example.com`).
3. Enter a secure Password – choose a strong password that you will remember.
4. Click the Sign Up button.

2.2 What Happens Next


- If your input meets the platform's requirements (unique email, valid format, and acceptable password), you will be automatically redirected to the login screen.
- If any of the fields contain invalid or duplicate information, the system will notify you to make corrections before you can proceed.

After successful registration, you can log in and begin using the platform to analyze chest X-rays.

3 Login

3.1 How to Log In

1. Navigate to the login page of the platform.
2. Enter your registered Email Address in the email field.
3. Enter your Password in the password field.
4. Click the Login button to access your account.



The image shows a login form for 'Scanalyze Imaging'. At the top is the logo, which consists of a blue square with a white stylized 'S' and the text 'Scanalyze Imaging' in blue. Below the logo are two input fields. The first is labeled 'Email' and contains the text 'steven@gmail.com'. The second is labeled 'Password' and contains six dots. Below these fields is a large blue button with the text 'Sign In' and a small cursor icon. At the bottom of the form is a link that says 'Don't have an account? Sign up'.

3.2 What Happens Next

- If the credentials are correct, you will be redirected to the Model Selection page, where you can choose the diagnostic tool you want to use.
- If the credentials are incorrect, an error message will appear stating, "Invalid email or password. Please try again."

3.3 Important Notes

- Ensure that you use the same email address and password you registered with.
- For security reasons, the platform enforces a maximum of 5 failed login attempts before temporarily locking the account. If this happens, wait 15 minutes before trying again or contact support for assistance.

4 Select Model

After logging in, you will be directed to the platform's Model Selection page. This page allows you to choose from several AI-powered diagnostic tools tailored for different types of medical imaging.

4.1 Available Models

The available models include:

- Chest X-Ray Analysis: Detects various chest conditions, including pneumonia, COVID-19, and other abnormalities.
- Brain MRI Segmentation: Identifies and segments brain tumors and lesions from MRI scans.
- Bone Fracture Detection: Detects fractures and abnormalities in X-ray images of bones.

4.2 How to Proceed

To analyze a chest X-ray image:

1. Locate the Chest X-Ray Analysis card.
2. Click on the "Chest X-Ray Analysis" option.

This action will take you to the CXR upload and analysis interface, where you can begin working with patient chest X-rays.

Medical Imaging AI Platform

Advanced AI-powered analysis for medical imaging, providing rapid and accurate diagnostics support



Chest X-Ray Analysis

Advanced AI model for detecting various chest conditions including COVID-19, pneumonia, and other abnormalities



Brain MRI Segmentation

Precise segmentation of brain MRI scans for tumor detection and analysis



Bone Fracture Detection

Rapid identification of bone fractures and abnormalities in X-ray images

5 Upload CXR Image

After selecting the Chest X-Ray Analysis model, you will be directed to the image upload page. This is where you submit the patient's chest X-ray (CXR) for automated analysis.

5.1 How to Upload an Image

The upload interface includes a clear image upload box, as shown below:

You can upload a medical image in either of the following ways:


- Drag and drop your CXR image file into the upload area.
- Or, click the upload box to manually browse and select an image file from your device.

Once an image is uploaded, the platform will begin processing the file using its AI model.

5.2 Important Notes

- Ensure the image is in a supported format (e.g., JPG, PNG, or DICOM).
- The AI system will generate predictions, a probability distribution of possible chest diseases, a detailed report of findings, and a visual explanation (e.g., heatmap) for interpretability.
- This tool is designed to assist with medical diagnostics and should not be used as a standalone diagnostic decision-maker.

A disclaimer is included at the bottom of the page stating that the tool is educational and should not replace professional clinical judgment.



Upload Medical Image
Drag and drop your image here, or click to select

00000211_036.png

Disclaimer: This tool is for educational purposes only and is not intended to replace professional medical diagnosis. This tool aims to assist such medical diagnosis, however, it should not be used in isolation.

For more information, visit:

- Radiopaedia
- CDC

6 Predicted Results and Report of Findings

After a chest X-ray image is uploaded, the platform automatically performs an AI-based analysis and generates a comprehensive diagnostic report. This report is intended to assist radiologists in interpreting the model's findings with full transparency and interactivity.

6.1 Report Overview

Upon completion of the analysis, a completed report is displayed, including the following key components:

- Timestamp indicating when the analysis was performed.
- A list of predicted findings based on the AI's highest confidence scores.
- Detailed medical descriptions of each predicted condition.
- A confidence slider to filter results.
- A probability distribution table of all considered diseases.
- A color-coded heatmap for quick visual interpretation of probabilities.

6.2 Confidence Slider

Located near the top of the report is the Confidence Threshold Slider, allowing users to adjust the minimum probability for a condition to be displayed.

- Sliding right filters out lower-confidence predictions and shows only those with high certainty.
- This feature allows radiologists to tailor their level of diagnostic confidence based on their clinical context, focusing only on the most probable conditions if desired.

6.3 Findings

Immediately below the confidence slider, the model's most likely conditions are listed under "Findings." These represent the diagnoses with the highest probabilities. For example, in the displayed report, the model has identified Effusion and Edema as the top findings.

Description of Predicted Conditions

Each predicted condition is accompanied by a detailed explanation, which includes:

- Clinical definition: Explaining what the disease is and its implications (e.g., edema is fluid accumulation often seen in congestive heart failure).
- Radiologic description: Highlighting how the disease typically appears on a chest X-ray (e.g., bathing opacities, Kerley B lines, or blunting of the costophrenic angle).

These descriptions help radiologists compare the AI findings with their own image interpretation, strengthening the diagnostic decision-making process.

6.4 Probability Distribution Table

The system provides a detailed probability distribution of all analyzed disease classes. Each condition is shown with its associated percentage, representing the model's confidence.

- This table is sorted to emphasize high-confidence predictions while also revealing lower-likelihood conditions for full transparency.
- It is particularly helpful when considering differential diagnoses and diagnosing conditions with low probability.

6.5 Heatmap Visualization

Beneath the probability table, a heatmap view provides a fast and intuitive way to interpret disease probabilities.

- Conditions are presented as color-coded boxes, ranging from red (high probability) to green (low probability).
- This allows for quick visual assessment of the disease, with higher probability conditions (e.g., Effusion at 66.5%) shown in red, while less likely ones (e.g., Emphysema at 7.85%) appear in green.
- The heatmap components this numerical data, offering an auto-glance summary that enhances situational awareness.

Chest X-Ray Analysis Report

3/28/2025, 8:34:19 PM

✓ Analysis Complete

Confidence Threshold: 0%



Findings: Edema, Effusion

Edema:

Pulmonary edema is fluid accumulation in the alveolar spaces and interstitium of the lungs, often due to congestive heart failure. On X-ray, you might see Kerley B lines, perihilar 'batwing' opacities, or diffuse haziness in severe cases.

Effusion:


A pleural effusion is fluid in the pleural space. On an X-ray, it can manifest as blunting of the costophrenic angle, a meniscus sign, or layering on a lateral decubitus view.

Class Probabilities (≥ 0%):

| | | | |
|----------------|--------|---------------------|--------|
| Atelectasis: | 39.91% | Cardiomegaly: | 36.08% |
| Consolidation: | 33.32% | Edema: | 64.98% |
| Effusion: | 65.50% | Emphysema: | 7.26% |
| Fibrosis: | 14.04% | Hernia: | 8.88% |
| Infiltration: | 39.11% | Mass: | 21.41% |
| Nodule: | 11.32% | Pleural_Thickening: | 31.55% |
| Pneumonia: | 28.07% | Pneumothorax: | 6.20% |

Heatmap Visualization

Low → High



7 Exploring the CXR Image

Once the AI has completed its analysis and generated the report, the original chest X-ray (CXR) image is displayed on the screen. This image viewer is equipped with interactive tools to allow detailed inspection by radiologists.

7.1 Interactive Tools for Image Exploration

Radiologists can explore the CXR image using the following built-in tools:

- **Zoom Functionality:** You can zoom in on specific areas of the X-ray for a closer inspection. This is especially helpful when reviewing image findings mentioned in the report, such as observing Kerley B lines, fluid levels, or lung opacities.
- **Drag-to-Pan:** After zooming in, you can click and drag the image to move across different regions of the scan. This allows you to focus on particular anatomical structures or abnormalities in detail, without losing your place.
- **Responsive Exploration:** The viewer responds smoothly to mouse scroll or trackpad pinch gestures for zooming, and click-drag gestures for panning, mimicking a diagnostic workstation experience.

7.2 Practical Use


These interactive tools are designed to support and enhance the interpretation of the report findings. For example:

- If the report highlights effusion, the radiologist can zoom in on the lower lung fields or costophrenic angles to confirm fluid levels.
- If edema is noted, the upper and mid-lung zones can be inspected more closely for classic patterns such as vascular redistribution or perihilar haze.

The goal is to give radiologists the control to visually verify and explore AI-generated predictions within the image itself.

Zoom & Pan Instructions

You can scroll or pinch to zoom into specific areas of the chest X-ray, and click-drag to pan around. This allows for a closer inspection of any region.



8 Download Report

After reviewing the full analysis, radiologists have the option to save a complete copy of the findings using the Download Report button.

8.1 Purpose of the Report

The downloadable report serves as a comprehensive summary of the AI analysis and is suitable for:

- Storing in patient records for documentation or follow-up.
- Sharing with colleagues for consultation or second opinion.
- Archiving for research or institutional case tracking.

8.2 What the Report Includes

When the Download Report button is clicked, a formatted PDF file is generated containing all key diagnostic information:

- Timestamp of when the analysis was completed.
- Predicted findings selected by the AI model.
- Detailed medical descriptions of the predicted conditions, including radiographic cues.
- The full probability distribution table showing the AI's confidence scores for all possible conditions.
- The confidence threshold setting used at the time of analysis.
- The heatmap visualization summarizing likelihoods in a color-coded format.

This PDF mirrors the information seen on-screen, ensuring nothing is lost or omitted in the export.

8.3 Using the Button

1. Once the analysis report is displayed, locate the Download Report button near the bottom of the interface.
2. Click the button to begin the download.
3. The report will be saved to your local system in PDF format, ready for use in clinical documentation or further review.

This feature ensures radiologists can retain and share a consistent, clear, and complete copy of the AI findings at any time.

You can download a PDF version of this analysis for your records.

Download PDF

9 Frequently Asked Questions (FAQ)

9.1 General Questions

9.1.1 What is the purpose of the Scanalyze Imaging Platform?

The platform is designed to assist radiologists by providing AI-powered diagnostic support for chest X-rays, including disease predictions, probability distributions, and visual explanations.

9.1.2 Is this platform a replacement for a radiologist's diagnosis?

No, the platform is a decision-support tool meant to complement professional clinical judgment, not replace it. Radiologists should always verify AI-generated predictions with their own experience and knowledge.

9.2 Usage Questions

9.2.1 What image formats are supported for upload?

The platform supports JPG, JPEG, and PNG formats.

9.2.2 How long does it take to analyze an image?

The analysis typically takes less than 5 seconds, depending on the image size and server load.

9.2.3 Can I use the platform on mobile devices?

Yes, the platform is optimized for both desktop and mobile browsers.

9.3 Technical Questions

9.3.1 What happens if the platform cannot process my image?

If the image is unsupported or corrupted, an error message will appear. Ensure the image is in a supported format and try again.

9.3.2 How secure is my data?

The platform uses AES-256 encryption for data storage and transmission, ensuring your data is secure.

9.3.3 Can I download the analysis results?

Yes, you can download a comprehensive PDF report of the analysis using the "Download Report" button.

9.4 AI and Diagnostic Questions

9.4.1 How accurate are the AI predictions?

The AI model has been trained on a large dataset (NIH chest xray dataset) which has over 100,000+ images. It has 70%+ AUC ROC score which means it can discriminate between an existence of a disease and the non-existence of one well. However, predictions should always be verified by a radiologist.

9.4.2 What diseases can the AI detect?

The AI can detect 13 different conditions: Atelectasis, Cardiomegaly, Consolidation, Edema, Effusion, Emphysema, Fibrosis, Infiltration, Mass, Nodule, Pleural Thickening, Pneumonia, and Pneumothorax.

9.4.3 What should I do if the AI gives low-confidence predictions?

If the confidence scores are low, consult a radiologist for further evaluation.